



Guldmann™

## Guldmann GL5.2 Lifter – Quiz Questions

**1. When is using the GL5.2 lift appropriate?**

- A. When a powered mobile working tool is indicated
- B. Recover patients from the floor after a fall
- C. Moving to and from bed, toilet, commode chair, or bedside chair
- D. When the floor and environment is even and level
- E. All of the above**

**Answer E**

**2. Which of the following are accurate features GL5.2?**

- A. Service light indicator to alert users to service needed
- B. Integrated weighing scale
- C. Battery can be charged on board by plugging in the GL5.2
- D. A and C**

**Answer D**

**3. What is the SWL (safe working load) of the GL5.2?**

- A. 155 kg. / 341 lbs.
- B. 205 kg. / 450 lbs.**
- C. 254 kg. / 560 lbs.
- D. 175 kg. / 385 lbs.

**Answer: B**

**4. When do you lock the GL5.2?**

- A. When lifting and lowering the patient
- B. When patients weight exceeds 275lbs
- C. When the GL5.2 is parked and not in use**
- D. All of the above

**Answer: C**

**5. The GL5.2 can be operated when plugged in to charge**

- A. true
- B. false**

**Answer: B**



Guldmann™

## Guldmann GL5.2 Active Lifter –

### Quiz Questions

6. Which are accurate statements regarding the GL5.2?

- A. GL5.2 can be used to transport patients long distances
- B. GL5.2 can be used for in bed care tasks such as rolling, repositioning, lift limbs
- C. GL5.2 has both low and high lifting capability
- D. B and C

Answer: D

7. Which item is NOT a responsible caregiver use of the GL5.2?

- A. Conducting a full risk assessment of user's functional communication and capability before use
- B. Explain to the user the use of and the procedural steps for using the GL.1 in advance of and during transfers and movement
- C. Leaving the user unattended when GL5.2 is in use
- D. Referring to the Product Manual and Manufacturer's literature for guidance

Answer: C

8. It is possible to tell what percentage charge the GL5.2 battery has.

- A. True
- B. False

Answer: A

9. What is maximum recharging time for GL5.2?

- A. 4 hours
- B. 2 hours
- C. overnight
- D. 3 hours

Answer: A

10. Name two ways to perform emergency lowering

- A. Push downward on the lifting boom arm of GL5.2
- B. Push the arrow down button on GL5.2's control box beneath the red emergency stop
- C. Manually rotate the red cylinder beneath the lifting boom arm
- D. Activate the base legs to their widest position
- E. B and C
- F. A and D

Answer: E