

## LP50 & LP50H Technical specifications

### Functions

Lifting capacity : up to 375 kg  
Lifting speed : 13 mm/sec.

### Operation

Electric : Operating panel/wireless operation  
Emergency stop : Electronic  
Safety cutout : Mechanical/electrical safety cutout

### Materials

Materials : Aluminium/galvanized steel  
Safety edge : Anodized aluminium with rubber strip  
Operating panel : Powder coated steel. Front rustproof steel.  
Platform floor finish : Optional. The weight of the floor finish will be deducted from the lifting capacity.

### Measurements

Thickness of floor finish : 22 mm  
Lifting height : 830 mm. Be aware of local regulations in relation to lifting height.  
Pit depth : 155 mm. Frame for embedment must always be supplied by Guldmann.  
Cantilever LP50H : 900 mm

### Electrical parts

On/off : Automatically  
Power consumption  
Stand-by : 3W  
Operation : 14.7 W  
Periodic load : 10/90  
Fuse : T10A  
Batteries : 2 x 12V/7.0 Ah  
Gas-proof and maintenance-free  
Power supply  
Primary : 230 V/115 V  
Secondary : 24 V DC

### Accessories/

**Miscellaneous** : External operating panel  
Wireless operation

### Marking

The product is manufactured in compliance with the Council Directive 2006/42/EC of May 17th 2006.

### Item Number/Model

13013-xxx : LP50 vertical lifting platform  
13012-xxx : LP50H cantilever vertical lifting platform

### Classification

According to ISO 9999 : 18 30 06

