



Guldmann Inc.
14401 McCormick Drive, Unit A
Tampa, FL 33626
Tel. 800 664 8834
Tel. 813 880 0619
Fax 813 880 9558
info@guldmann.net
www.guldmann.net

GH Cross hanger **Technical specifications**

Functions

Lifting capacity, max. : 1100 lbs (500 kg)
Suspension points : 4

Weight and Materials

Weight, totally : 19.4 lbs (8.8 kg)
Lifting hanger : Rustproof steel, powder coated

Measurements

A : 33³/₄" (860 mm)
B : 21¹/₄" (540 mm)
C : 10" (250 mm)

Accessories/Miscellaneous



Lifting sling : See brochure

Marking

The product is manufactured in compliance with the Council Directive 93/42/EEC of June 14th 1993, including amendments, as medical device class 1.

Classification

according to ISO 9999 : 12 36 03



| Time to care |

Guldmann™