



Guldmann Inc.
Tel. 800 664 8834
Tel. 813 880 0619
info@guldmann.net
www.guldmann.net

H-Hanger Technical specifications

Functions

Lifting capacity, max. : 825 lbs (375 kg)
Suspension points : 2

Weight and Materials

H-Hanger : 3.5 lbs (1,6 kg)
Lifting hanger : Rustproof steel covered
with polyurethane foam
Suspension points : Rustproof steel

Labelling

The product is manufactured in compliance with regulation (EU) 2017/745 of the European parliament and of the Council of 5 April 2017, as medical device Class I.

Item Number/Model

556950 : H-Hanger

Classification



according to ISO 9999 : 12 36 03

Accessories/Miscellaneous

Swing kit : 554903
Lifting hanger
Medium : 556890
Small : 556880
X-Small : 556970

H-Hanger + Lifting hanger

556950 + 556890 : C = 13 1/2" (340 mm)
556950 + 556880 : C = 12 1/4" (312 mm)
556950 + 556970 : C = 11 1/2" (289 mm)



Time to care

Guldmann™