



## Lifting hanger (2016) Technical specifications

### Functions

Lifting capacity, max. : 375 kg  
Suspension points : 4

### Weight and Materials

Medium, totally : 1,8 kg  
small, total : 1,6 kg  
X-small, total : 1,35 kg

Lifting hanger : Rustproof steel covered  
with polyurethane foam  
Suspension points : Rustproof steel

### Measurements

556890 medium

A : 600 mm  
B : 534 mm  
C : 128 mm

556880 small

A : 500 mm  
B : 440 mm  
C : 100 mm

556870 x-small

A : 395 mm  
B : 343 mm  
C : 77 mm

### Accessories/Miscellaneous


Lifting sling : See brochure

### Labelling

The product is manufactured in compliance with regulation (EU) 2017/745 of the European parliament and of the Council of 5 April 2017, as medical device Class I.

### Classification

according to ISO 9999 : 12 36 03



Time to care