



**Guldmann Inc.**  
Tel. 800 664 8834  
Tel. 813 880 0619  
Fax 813 880 9558  
info@guldmann.net  
www.guldmann.net

## Lifting hanger (2016) Technical specifications

### Functions

Lifting capacity, max. : 825 lbs (375 kg)  
Suspension points : 4

### Weight and Materials

Medium, totally : 4 lbs (1.8 kg)  
SMALL, total : 3.5 lbs (1.6 kg)  
X-SMALL, total : 3 lbs (1.35 kg)

Lifting hanger : Rustproof steel covered  
with polyurethane foam  
Suspension points : Rustproof steel

### Measurements

556890 medium

A : 24" (600 mm)  
B : 21 3/8" (534 mm)  
C : 5 1/8" (128 mm)

556880 small

A : 20" (500 mm)  
B : 17 1/2" (440 mm)  
C : 4" (100 mm)

556870 x-small

A : 15 3/4" (395 mm)  
B : 13 3/4" (343 mm)  
C : 3" (77 mm)

### Accessories/Miscellaneous


Lifting sling : See brochure

### Labelling

The product is manufactured in compliance with regulation (EU) 2017/745 of the European parliament and of the Council of 5 April 2017, as medical device Class I.

### Classification

according to ISO 9999 : 12 36 03



Time to care

Guldmann™