



V. Guldmann A/S
Tel. +45 8741 3100
Fax +45 8741 3131
info@guldmann.com
www.guldmann.com

Free-standing room covering rail system, height adjustable

Technical specifications

Functions

Max. capacity : 350 kg

Weight and Materials

Free-standing room covering

rail system, 3 x 3 m : 151 kg
4 Vertical poles : 11 kg, each
4 B-rails : 16,4 kg, each
4 top corner brackets : 4,8 kg, each
1 Traverse rail
with trolleys : 22,2 kg

Free-standing room covering

rail system, 4 x 4 m : 178 kg
4 Vertical poles : 11 kg, each
4 B-rails : 21,8 kg, each
4 top corner brackets : 4,8 kg, each
1 Traverse rail
with trolleys : 27,6 kg

Rails and poles : Anodized aluminium
Brackets : Powder varnished steel

Measurements

Height min. : 2040 mm
Height max. : 2840 mm

Length, 3 x 3 m : 3200 mm
Width, 3 x 3 m : 3200 mm




Length, 4 x 4 m : 4200 mm
Width, 4 x 4 m : 4200 mm

Labelling

The product is manufactured in compliance with the Council Directive 93/42/EEC of June 14th 1993, including amendments, as medical device class 1.

Item Number/Model

555262 : Free-standing room covering
rail system, height adjustable - 3 x 3 m
555570 : Free-standing room covering
rail system, height adjustable - 4 x 4 m



Time to care