



Guldmann Inc.
 Tel. 800 664 8834
 Tel. 813 880 0619
 Fax 813 880 9558
 info@guldmann.net
 www.guldmann.net

Free-standing room covering rail system, height adjustable

Technical specifications

Functions

Max. capacity : 770 lbs (350 kg)

Weight and Materials

Free-standing room covering

rail system, 3 x 3 m : 332 lbs (151 kg)
 4 Vertical poles : 24 lbs (11 kg), each
 4 B-rails : 36 lbs (16,4 kg), each
 4 top corner brackets : 11 lbs (4,8 kg), each
 1 Traverse rail
 with trolleys : 49 lbs (22,2 kg)

Free-standing room covering

rail system, 4 x 4 m : 392 lbs (178 kg)
 4 Vertical poles : 24 lbs (11 kg), each
 4 B-rails : 48 lbs (21,8 kg), each
 4 top corner brackets : 11 lbs (4,8 kg), each
 1 Traverse rail
 with trolleys : 61 lbs (27,6 kg)

Rails and poles : Anodized aluminium
 Brackets : Powder varnished steel

Measurements

Height min. : 2040 mm
 Height max. : 2840 mm

Length, 3 x 3 m : 3200 mm
 Width, 3 x 3 m : 3200 mm




Length, 4 x 4 m : 4200 mm
 Width, 4 x 4 m : 4200 mm

Labelling

The product is manufactured in compliance with regulation (EU) 2017/745 of the European parliament and of the Council of 5 April 2017, as medical device Class I.

Item Number/Model

555262 : Free-standing room covering
 rail system, height adjustable - 3 x 3 m
 555570 : Free-standing room covering
 rail system, height adjustable - 4 x 4 m



Time to care