



## Positioning Lock Technical specifications

### Function

Usage : To lock GH3 hoist and traverse rail  
in a certain position

### Weight

Lock : 3 x 1,4 kg

### Materials

Lock function : Rustproof and galvanized steel  
Covers : Impact-resistant UL 94 V-0 flame  
retardant recyclable plastic

### Electronics

On/off : Automatic when used.  
LED indicator : Operation indicator  
Power supply : 33 V AC, 2.0 A  
Supply voltage,  
transformer : 115/230 V AC  
Battery : 1 x 24V NiMH  
Max charge time : 2 hours

### Measurements


A : 120 mm  
B : 150 mm  
C : 70 mm  
D : min. 10 mm clearance

### Markings

The product is manufactured in compliance with regulation (EU) 2017/745 of the European parliament and of the Council of 5 April 2017, as medical device Class I.

### Item Number/Model

557600 : Positioning Lock



| Time to care |