



Guldmann Inc.
Tel. 800 664 8834
Tel. 813 880 0619
Fax 813 880 9558
info@guldmann.net
www.guldmann.net

Positioning Lock Technical specifications

Function

Usage : To lock GH3 hoist and traverse rail
in a certain position

Weight

Lock : 3 x 3 lb (1,4 kg)

Materials

Lock function : Rustproof and galvanized steel
Covers : Impact-resistant UL 94 V-0 flame
retardant recyclable plastic

Electronics

On/off : Automatic when used.
LED indicator : Operation indicator
Power supply : 33 V AC, 2.0 A
Supply voltage,
transformer : 115/230 V AC
Battery : 1 x 24V NiMH
Max charge time : 2 hours

Measurements


A : 4 3/4" (120 mm)
B : 6 (150 mm)
C : 2 3/4" (70 mm)
D : min. 1/2" (10 mm) clearance

Markings

The product is manufactured in compliance with regulation (EU) 2017/745 of the European parliament and of the Council of 5 April 2017, as medical device Class I.

Item Number/Model

557600 : Positioning Lock



Time to care