



Multi-flexible design solutions for residential nursing homes

There is great focus today on building flexible residential nursing homes that allow us to meet the needs and wishes of each resident while also creating the best working conditions for caregivers.

The needs of individual residents vary greatly and can change over time, e.g. from being ambulatory to a more reliant and possibly

bedridden nature. Furthermore, requirements for accommodation can change when a new resident moves in.

Guldmann and Flexi Foldedøre's solution of mobile walls and revolving/folding doors with ceiling hoist system allows for a multi-flexible design of residential nursing homes with additional options for

changing the design when new needs arise, without dismantling or rebuilding.

It's about using the best building solutions to create the best possible conditions for a growing number of older residents and not least caregivers.



1. Combining Flexi's revolving/folding doors and mobile walls with a Guldmann GH3 patient lift system provides unhindered access to the entire residence. Parallel rails are built discreetly into the ceiling to function as tracks for a room covering system. This solution provides a flexible division and design of the residential space.

2. Flexi's Cubitus revolving/folding doors provide the user with easier access to the bathroom, for example, since the doors provide maximum free width and can easily be operated from a wheelchair.

3. Flexi's mobile walls allow the resident to divide the rooms of the residence into different sizes according to preference and needs.

4. Folding doors and mobile walls run effortlessly into transverse Guldmann ceiling hoist rails. Both folding doors and mobile walls are operated using a transverse rail crossing module out to the side so that a GH3 ceiling hoist can cover and serve the entire residence unhindered.

Guldmann™

V. Guldmann A/S
Graham Bells Vej 21-23A
DK-8200 Aarhus N
Tel. +45 8741 3100
info@guldmann.com
www.guldmann.com



FLEXI Foldedøre ApS
Skt. Georgs Vej 1
DK-8260 Viby J.
Tel. +45 8628 4066
flexi@flexi-foldedoere.dk
www.flexi-foldedoere.dk