



Guldmann Inc.
Tel. 800 664 8834
Tel. 813 880 0619
Fax 813 880 9558
info@guldmann.net
www.guldmann.net

Horizontal lifting support with weight adjustment, 825 lbs Technical specifications

Functions

Lifting capacity, max. : 825 lbs
Suspension points : 8

Weight and Materials

Total weight : 25,5 lbs
Lifting hanger : Powder coated steel, RAL 9016

Measurements

A : 42 inch
B : 4 ½ inch
C : 55 ¾ inch
D : 20 ¼ inch
E : 19 ¾ inch
F : 21 inch
G : 22 ¾ inch
H : 28 ½ inch
I : 8 ½ inch

Accessories/ Miscellaneous

Lifting sling : See brochure

Labelling


The product is manufactured in compliance with regulation (EU) 2017/745 of the European parliament and of the Council of 5 April 2017, as medical device Class I.

Item Number/Model

28456 : Horizontal lifting support with weight adjustment

Classification

according to ISO 9999 : 12 36 03



Time to care

Guldmann™