



Guldmann Inc.
 Tel. 800 664 8834
 Tel. 813 880 0619
 Fax 813 880 9558
 info@guldmann.net
 www.guldmann.net

Free-standing room covering rail system, adjustable

Technical specifications

Functions

Max. capacity : 560 lbs (255 kg)

Weight and Materials

Free-standing room covering

rail system, 3 x 3-4 m : 352 lbs (160 kg)
 4 Vertical poles : 24 lbs (11 kg), each
 4 B-rails : 36 lbs (16,4 kg), each
 4 top corner brackets : 15 lbs (7 kg), each
 1 Traverse rail
 with trolleys : 49 lbs (22,2 kg)

Free-standing room covering

rail system, 4 x 4-5 m : 411 lbs (187 kg)
 4 Vertical poles : 24 lbs (11 kg), each
 4 B-rails : 48 lbs (21,8 kg), each
 4 top corner brackets : 15 lbs (7 kg), each
 1 Traverse rail
 with trolleys : 61 lbs (27,6 kg)

Rails and poles : Anodized aluminium
 Brackets : Powder varnished steel

Measurements

Height min. : 2040 mm
 Height max. : 2840 mm

Free-standing room covering rail system, 3 x 3-4 m

Width : 3200 mm
 Length : 3200-4200 mm
 Working area : 3000 x 3000 mm

Free-standing room covering rail system, 4 x 4-5 m




Width : 4200 mm
 Length : 4200-5200 mm
 Working area : 4000 x 4000 mm

Labelling

The product is manufactured in compliance with regulation (EU) 2017/745 of the European parliament and of the Council of 5 April 2017, as medical device Class I.

Item Number/Model

554729 : Free-standing room covering rail system, adjustable - 3 x 3-4 m
 557968 : Free-standing room covering rail system, adjustable - 4 x 4-5 m



Time to care