



Free-standing room covering rail system, adjustable

Technical specifications

Functions

Max. capacity : 255 kg

Weight and Materials

Free-standing room covering

rail system, 3 x 3-4 m : 160 kg
4 Vertical poles : 11 kg, each
4 B-rails : 16,4 kg, each
4 top corner brackets : 7 kg, each
1 Traverse rail
with trolleys : 22,2 kg

Free-standing room covering

rail system, 4 x 4-5 m : 187 kg
4 Vertical poles : 11 kg, each
4 B-rails : 21,8 kg, each
4 top corner brackets : 7 kg, each
1 Traverse rail
with trolleys : 27,6 kg

Rails and poles : Anodized aluminium
Brackets : Powder varnished steel

Measurements

Height min. : 2040 mm
Height max. : 2840 mm

Free-standing room covering

rail system, 3 x 3-4 m
Width : 3200 mm
Length : 3200-4200 mm
Working area : 3000 x 3000 mm

Free-standing room covering




rail system, 4 x 4-5 m
Width : 4200 mm
Length : 4200-5200 mm
Working area : 4000 x 4000 mm

Labelling

The product is manufactured in compliance with regulation (EU) 2017/745 of the European parliament and of the Council of 5 April 2017, as medical device Class I.

Item Number/Model

554729 : Free-standing room covering
rail system, adjustable - 3 x 3-4 m
557968 : Free-standing room covering
rail system, adjustable - 4 x 4-5 m



Time to care