



**V. Guldmann A/S**  
 Tel. +45 8741 3100  
 Fax +45 8741 3131  
 info@guldmann.com  
 www.guldmann.com

## Free-standing room covering rail system, adjustable

### Technical specifications

#### Functions

Max. capacity : 255 kg

#### Weight and Materials

*Free-standing room covering*

*rail system, 3 x 3-4 m* : 160 kg  
 4 Vertical poles : 11 kg, each  
 4 B-rails : 16,4 kg, each  
 4 top corner brackets : 7 kg, each  
 1 Traverse rail  
 with trolleys : 22,2 kg

*Free-standing room covering*

*rail system, 4 x 4-5 m* : 187 kg  
 4 Vertical poles : 11 kg, each  
 4 B-rails : 21,8 kg, each  
 4 top corner brackets : 7 kg, each  
 1 Traverse rail  
 with trolleys : 27,6 kg

Rails and poles : Anodized aluminium  
 Brackets : Powder varnished steel

#### Measurements

Height min. : 2040 mm  
 Height max. : 2840 mm

*Free-standing room covering*

*rail system, 3 x 3-4 m*  
 Width : 3200 mm  
 Length : 3200-4200 mm  
 Working area : 3000 x 3000 mm

*Free-standing room covering*




*rail system, 4 x 4-5 m*  
 Width : 4200 mm  
 Length : 4200-5200 mm  
 Working area : 4000 x 4000 mm

#### Labelling

The product is manufactured in compliance with the Council Directive 93/42/EEC of June 14th 1993, including amendments, as medical device class 1.

#### Item Number/Model

554729 : Free-standing room covering rail system, adjustable - 3 x 3-4 m  
 557968 : Free-standing room covering rail system, adjustable - 4 x 4-5 m



Time to care