

To whom it may concern

V. Guldmann A/S
Graham Bells Vej 21-23A
DK-8200 Aarhus N, Denmark
Tel. +45 8741 3100
Fax +45 8741 3131
info@guldmann.com
www.guldmann.com
CVR. nr. 27 70 67 46

November 29, 2017

Declaration recycling and use of recycled materials in Guldmann A/S products.

Rails:

Material: 100 % aluminium, consisting of 20- 25 % recycled aluminium.
The material is 100 % recyclable.

Cover for rails:

Material: 100 % virgin polypropylene.
The material is 100 % recyclable.

Brackets:

Material: 100% steel (alloy), consisting of 20 % recycled steel.
The material is 100 % recyclable.

Screws and bolts:

Material 1: 100 % aluminium, consisting of 20 – 25 % recycled aluminium.
2: 100 % steel (alloy), consisting of 20 % recycled steel.
The material is 100 % recyclable.

Cover for hoist:

Material: 100 % virgin ABS
The material is 100 % recyclable.

PCB:

Material (average weight % per kg PCB):

- Gold: 0,073 %, consisting of up to 98 % recycled Gold
- Palladium: 0,042 %, consisting of up to 98 % recycled Palladium
- Silver: 0,76 %, consisting of up to 90 % recycled Silver
- Tin: 1,02 %, consisting of up to 20 % recycled Tin
- Antimony: 0,23 %, not recycled
- Copper: 22,43 %, consisting of up to 98 % recycled copper
- Nickel: 0,18 %, not recycled
- Iron: 2,13 %, not recycled
- Aluminium: 0,14 %, consisting of up to 20 % recycled aluminium
- Plastics and ceramics are currently not recycled because the values of these are insignificant compared to above mentioned metals.

Electrical cords:

Material: 100 % copper with 30 – 50 % recycled copper.
The material is 100 % recyclable.

Slings:

Material: 100 % virgin polyester.
The material is 100 % recyclable.

Belts:

Material: 100 % virgin polyester.
The material is 100 % recyclable.

Yours sincerely

V. Guldmann A/S



Brian Hansen
Chief Operations Officer