



## Cross hanger, 400 kg Technical specifications

### Functions

Lifting capacity, max. : 400 kg  
Suspension points : 4

### Weight and Materials

Weight, totally : 4,6 kg

Lifting hanger : Rustproof steel, powder coated

### Measurements

A : 192 mm  
B : 309 mm  
C : 520 mm

### Accessories/Miscellaneous

Lifting sling : See brochure

### Marking


The product is manufactured in compliance with regulation (EU) 2017/745 of the European parliament and of the Council of 5 April 2017, as medical device Class I.

### Item Number/Model

561610 : Cross hanger, 400 kg

### Classification

according to ISO 9999 : 12 36 03



Time to care