



GLS5.2

Active lifter

– an easy, comfortable and safe way to raise and support patients with partial or complete control of upper body

One lift for many functions

- Sit to stand
- Standing balance (with foot tray removed)
- Transfers from chair to chair
- Transfers from bed to chair
- Transfers from bed to toilet

Adjusts easily to comfortably fit most heights and sizes

- Padded and adjustable knee support
- Multiple sling attachment points
- Ergonomically profiled push handle for all caregiver sizes
- Multiple sling sizes to create optimum comfort and support (option)
- Add-on hip sling to assist with hips and buttocks raising (option)

Easy and safe to use and understand

- Electric up, down, and leg width adjustment
- Dual operations options – hand control or control box
- Rear foot push bar to increase maneuverability
- Manually and electrically operated emergency lowering
- Easily accessible emergency stop



Strong, reliable and powerful

- SWL 155 or 205 kg
- Designed for intensive use in professional environments
- Meets and exceeds all current norms and regulations in industry
- Strong 25,2 V / 2,9 Ah Li-Ion battery ensures long usage time, even in heavy use environments
 - Shorter charging time to full capacity
 - Prolonged battery life
 - More lifts from a single charge
- Strong, noiseless and fast Linak actuators
- Varied charging options – on board lift and off board lift via extra battery
 - Built in and easily accessible docking station for battery
 - Extra battery packs as option
 - Fully automatic switch mode charger
- Removable and easy to clean foot tray
- The battery capacity is indicated in 3 intervals in the hand control;
 - Low 0 - 25%
 - Medium 25 - 50%
 - High 50 - 100%Audio signal at low battery voltage when approx. 25%.
- Visual service indicator lights up when maintenance is required.

| Time to care |

V. Guldmann A/S
Tel. +45 8741 3100
info@guldmann.com
www.guldmann.com